

OMi

Motion Activated Experiences For Sensory And Learning Environments

Origin: UK



Products At A Glance



Use and Benefits

- Create a calming space
- Broaden awareness & sense of self
- Enrich language and communication
- Improve coordination/motor skills
- Encourage exploration and discovery
- Reward with stunning audio visuals
- Experience collaborative learning
- Smart view - connect any compatible HDMI device for seamless viewing experiences

1. omiVista install - interactive floor projection
2. omiVista Mobii - interactive table/floor projection
3. omiReflex - interactive wall projection
4. omiBeam - interactive sensory lighting
5. omiSky - interactive ceiling projection
6. omiScent - aroma delivery system
7. omiLED - colour wash lighting
8. omiControl - room control system

omiVista Mobii

Interactive Floor / Table Projection


Origin: UK
Certificate: IEC 62368-1

The omiVista install and Mobii systems can be projected onto any horizontal surface. Such as table, tent, wheelchair tray, trampoline and there are many more creative ways to use the system including beds, bedside tables, flooring and pale carpets, etc.



Features

- Every omiVista system comes pre-programmed with the Sensory Fun Suite (Over 300 ready-made applications)
- The system is sensitive, responds to the smallest movement
- Editing and creation of new content made easy with the Window based user interface supporting most popular file formats for images, videos and sound
- Simple navigation of content using only 5 buttons and volume control



New Updates:

- A new and improved ergonomic Bluetooth remote control
- New improved movement tracking technology
- Easy access focus controls now located on top panel
- Laser Projector for clearer image and low maintenance

Enabling users to :

- Import their own images, videos and music
- Create person - centred Apps for even greater engagement



Care Suite

- Scenes and sounds from nature including therapeutic colouring
- Relaxing virtual water with responsive fish and the sounds of the sea
- Interactive quizzes to encourage/maintain memory and recall
- Physical games and activities to encourage participation/develop reflexes
- Nostalgic themes and music to stimulate discussion and reminiscence therapy
- Social games for shared enjoyment

Multi-Sensory

omiBudii

Sensory System for SEND

Origin: UK

By detecting the smallest of movements, Budii provides sensory stimulation, supports learning through fun, improves communication and motor skills, reduces anxiety, and has a calming effect. The projection technology turns any table or floor surface into an exciting and vibrant play area. The unit is easily wall or shelf mounted for in-home use.



Budii is specifically designed for home users with special education needs and disabilities (SEND). Budii supports the whole range of special needs in SEND, and has proven especially helpful to those with:

- Spina bifida
- Autism (ASD/ASC)
- Weak/poor muscle tone
- Cerebral palsy
- Balance and co-ordination difficulties
- Muscular dystrophy



Specifications

Dimensions (L x W x H)	13.5cm x 23cm x 26cm
Unit Weight	4.5kg
Recommended mounting height	1.4m - 2.0m (without the need to recalibrate camera)
Software	Non-editable (Containing 120 activities), 3 different standard suites: Budii Sensory, Budii ASN, Budii Care
Budii Flip Stand	Adjustable up to 2.1 meters from the floor



Budii Flip Stand

The Budii electric swivel mobile stand allows effortless height adjustment with foot-operated buttons. Featuring a swivel bracket, it enables projection onto tables, floors, walls, and ceilings.

Stand Height With Budii: from 1.1m to 2.1m fully extended

Image Size On Floor: up-to 2m x 1.2m wide

Complete Weight: 30kg approximately (including Budii)

Budii Manual Mobile Stand

The Budii manual mobile stand offers easy, precise height adjustment by turning the crank handle. Perfect for both table and floor projection.

Stand Height With Budii: from 1.4m to 2.2m fully extended

Image Size On Floor: up-to 2m x 1.2m wide

Complete Weight: 23kg approximately (including Budii)

omiReflex mobile

Interactive Wall Projection

Origin: UK



The omiReflex is a system that reveals a moving mirror image of the user with 'floating' objects that they can interact with. e.g. scattering leaves/petals, popping bubbles and bouncing spacemen.

It can project onto any pale vertical surface.

The omiReflex+ version is fully customisable allowing users to import their own images, sounds and videos.

Use and Benefits

- Coordination/motor skill development
- Proprioceptive and vestibular feedback
- Therapeutic calming activities
- Exercise and posture
- Social enjoyment and communication

Specifications

omiReflex Content Suite with 100+ apps

Wireless Remote Control

Approximate Projection Dimensions:
 0cm from the wall 1m x 0.75m,
 20cm from the wall 1.5m x 1.1m,
 35cm from the wall 2m x 1.5m




OMi Games

Origin: UK


OMi Games are a great enhancement for omiVista, omiReflex, and omiWall systems. Games can be played on table, floor, or wall surfaces.

Catch Games



- Basketball
- Bird Feeding Game
- Coconut Gorilla
- Golf Game

Ball Games



Scoring or Non-Scoring, your game, your way!

- American Football
- Basketball
- Basketball - Non-Scoring
- Beachball

Double Spat



Experience Catch Games in two thrilling ways – Play side by side or opposite each other.

- Fantasy
- Gem Mining
- Greek Mythology
- Nature

Sort Pro



Organise and categorise different items, from letters to colours, and challenge your mind in an engaging way.

- Farm Animals
- Food and Countries
- Rooms of the House
- Sort Capital Letters

Puzzle



Challenge yourself to complete colourful puzzles that sharpen focus and enhance cognitive abilities

- Nature
- Landmarks
- Biomes
- Cities

Match Games



Test your memory and matching skills with vibrant tiles that engage and challenge your mind.

- Beach
- Fruit
- Silhouettes
- Transport

Quiz Pro Hong Kong Version NEW

Drawing inspiration from the city's love for trivia, Quiz Pro HK Version features a localized question bank of over 500 unique challenges to capture the pulse of Hong Kong, including:

- Local Sports & Icons
- Canto-Pop & Cinema
- History & Heritage
- Hong Kong General Knowledge


Vecta Delight

Mobile Sensory System

Origin: UK



Vecta - a portable, versatile sensory experience for any space. The Vecta Distraction Station is ideal where space is at a premium, such as in hospitals and nursing homes. Converts waiting rooms, treatment rooms and bedrooms into relaxing, distracting and empowering Multi-Sensory Environments. The stable Vecta Distraction Station can be moved easily from room to room on its sturdy wheels and requires only one plug to transform the room with the 'Big 3' Multi-Sensory items: Bubble Column, Projector and Fibre Optics.



Product code: 9VCTAH

Dimensions: 150cm x 45cm x 58cm